

# A Web-Based Intelligent System for Automated Detection, Classification, and Analysis of Microplastics from Microscopic Images Using Ultralytics YOLO

Perarasan K<sup>1</sup>, Dr. K. Santhi<sup>2</sup>

Department Of Information Technology, Dr. N.G.P Arts and Science College, Coimbatore, Tamil Nadu, India<sup>1</sup>

Professor, Department of Information Technology, Dr. N.G.P Arts and Science College, Coimbatore, Tamil Nadu, India<sup>2</sup>

**Abstract:** Microplastics (MPs (Micro, (Microplastics), 1  $\mu\text{m}$ –5 mm) and nano plastics (NPs(Nanoplastics), <1  $\mu\text{m}$ ) pose significant environmental and health risks, yet traditional detection methods are slow and resource-intensive. MPWebAI is a browser-based platform enabling users to upload microscopic images for automated MP detection, morphological classification (fibre, film, fragment, pellet), and quantitative analysis using Ultralytics YOLOv11m. Trained on a 7,200-image curated dataset, the model achieves mAP@50 of 95.4%, precision 94.2%, and recall 92.8%. The Fast API + React web app provides annotated images, particle counts, size statistics, and reports in seconds, requiring no specialized hardware beyond a standard microscope and smartphone. A review of current NP detection methods—including optical sieves, SERS, nano-DIHM, and AI-enhanced spectroscopy—highlights challenges for sub-micron particles such as matrix interference, spectral overlap, and the absence of certified reference materials. The paper outlines extensions integrating higher-resolution imaging and hybrid AI-spectroscopy. MPWebAI democratizes monitoring for researchers and citizen scientists in resource-limited regions. The system is open-source.

**Keywords:** Microplastics, Nano plastics, Ultralytics YOLOv11, Object Detection, Web-Based AI, Environmental Monitoring, Deep Learning.

## I. INTRODUCTION

Plastic pollution has reached critical levels, with microplastics (MPs) and the even smaller nano plastics (NPs, typically <1  $\mu\text{m}$ , often <100 nm) now detected in virtually every environmental compartment and human tissue. NPs are particularly concerning due to their ability to cross biological barriers, high surface-area-to-volume ratio, and potential for chemical leaching. Recent studies have detected NPs in human blood, placenta, and breast milk, raising profound concerns about their toxicological effects, including oxidative stress, immune disruption, and potential associations with human illnesses.

While MP detection has benefited from deep learning, NP detection remains challenging because of size, low concentrations, and matrix interference. Conventional analytical techniques such as Fourier-transform infrared spectroscopy (FTIR), Raman spectroscopy, and pyrolysis-gas chromatography-mass spectrometry (Py-GC-MS) are limited by high cost, labour-intensive sample preparation, and poor sensitivity at the nanoscale. The diffraction limitation of FTIR makes it intrinsically challenging to detect particles smaller than 1 micrometre, positioning Raman spectroscopy as a more promising but still methodologically difficult alternative.

MPWebAI provides an immediate, accessible solution for MPs and lays the foundation for NP analysis by supporting higher-magnification uploads and planning multi-modal extensions. The platform addresses the urgent need for scalable monitoring solutions that can operate beyond advanced laboratories, where detection of sub-micron particles has remained nearly impossible.

### Key Contributions

First open, web-deployed YOLOv11-based multi-class MP morphology detector with integrated sizing.

Comprehensive review of state-of-the-art NP detection methods (2024–2026), including optical sieves, SERS substrates, and AI-enhanced spectroscopy.

Roadmap for extending the platform to nano plastics using emerging optical, spectroscopic, and AI techniques.

**II. LITERATURE REVIEW**

Detection of microplastics (MPs, 1–5,000  $\mu\text{m}$ ) and nano plastics (NPs, 94% accuracy), electroanalytical sensors, microfluidic platforms, hyphenated Spectro-microscopic techniques, and mass spectrometry innovations like Hold-MS for direct airborne NP analysis. Persistent challenges involve the absence of certified NP reference materials, contamination risks, differentiation from natural colloids, and spectral overlaps, driving calls for multimodal, AI-integrated workflows. Overall, trends favor accessible, scalable tools merging imaging, spectroscopy, and nanotechnology to enable widespread monitoring beyond specialized laboratories.

**III. ROLE OF YOLO IN MICROPLASTIC DETECTION**

The Ultralytics YOLO (You Only Look Once) framework plays a central role in enabling fast and accurate automated detection of microplastics from microscopic images. Unlike traditional multi-stage detectors, YOLO performs object localization and classification in a single forward pass, making it highly suitable for real-time environmental analysis applications.

**A. Single-Stage Object Detection**

YOLO treats object detection as a regression problem and processes the entire image in one pass through the neural network. The input image is divided into a grid, and for each grid cell the model predicts bounding box coordinates, object confidence score, and class probabilities. This unified approach significantly reduces computational overhead compared to two-stage detectors and enables near real-time inference, essential for high-throughput microplastic analysis.

**B. Accurate Localization of Microplastics**

Microplastics in microscopic images often vary in size, shape, and orientation. YOLO's multi-scale feature extraction allows the model to detect thin fibres, irregular fragments, transparent films, and spherical pellets. The feature pyramid network within the YOLO architecture helps capture both small and large particles, improving detection robustness in complex backgrounds.

**C. Morphological Classification Capability**

Beyond detection, YOLO simultaneously classifies each detected object into predefined morphological categories: fibre, fragment, film, and pellet. This eliminates the need for separate classification pipelines and enables an end-to-end automated workflow.




Fig 1: Plastic detection

**D. Real-Time Processing Performance**

One of the major advantages of YOLO is its high inference speed. The lightweight architecture and optimized convolutional layers enable the system to process microscopic images within milliseconds. This supports rapid batch analysis, real-time laboratory workflows, scalable web deployment, and interactive user feedback—performance difficult to achieve using conventional image processing methods.

**E. Robustness to Complex Backgrounds**

Microscopic samples often contain noise, organic debris, and uneven illumination. YOLO's deep convolutional backbone learns discriminative features that help distinguish true microplastic particles from background artifacts. Data augmentation during training further improves the model's generalization ability.

**F. Integration with the Web-Based Pipeline**

The Ultralytics implementation provides a streamlined API that can be easily integrated with Python-based web frameworks. In the proposed system, the YOLO model is deployed on the server side, automatically processing uploaded images and returning detection results. This tight integration enables automated inference, scalable deployment, minimal user intervention, and efficient resource utilization.




Fig 2: YOLO Architecture

**G. Summary of Benefits**

The use of YOLO provides high detection accuracy, real-time inference speed, end-to-end detection and classification, robust performance on microscopic data, easy deployment in web environments, and scalability for large datasets.

**IV. EMERGING AI-INTEGRATED & HYBRID METHODS**

Recent advancements in NP detection demonstrate a clear trend toward hybrid methodologies that combine physical capture mechanisms, spectroscopic analysis, and machine learning. These innovations address the limitations of traditional methods and offer pathways for integration with platforms like MPWebAI.

**A. Optical Sieve Technology**

Researchers at the University of Stuttgart and the University of Melbourne have developed an "optical sieve"—a test strip with tiny cavities etched into a semiconductor substrate that captures nano plastic particles by size. When particles fall into Mie voids (200 nm to 1 μm), distinct colour changes occur that are visible under an ordinary optical microscope. This method enables rapid, low-cost detection without requiring electron microscopy or trained personnel. The technology has been validated using lake water samples spiked with nano plastics at concentrations of 150 μg/ml, demonstrating its potential for field deployment.

**B. AI-Assisted Nano-Digital In-line Holographic Microscopy (nano-DIHM)**

A cutting-edge nano-DIHM system developed for real-time and in-situ MNP research incorporates deep neural network functionality for rapid detection in environmental waters. This 4D (3D + time) tracking capability measures sedimentation velocities of plastics and distinguishes various polymer types (PE, PP, PS, PET, PVC, PUR) from organic and inorganic materials. The integration of deep learning enables automated classification and quantification, addressing the challenge of distinguishing NPs from natural colloids.

**C. Machine Learning-Enhanced Raman Spectroscopy**

Raman spectroscopy faces a fundamental trade-off between analysis speed and signal-to-noise ratio (SNR). Recent work demonstrates that machine learning can overcome this limitation: a Bi-CNN with attention mechanism achieved 99.29% accuracy in classifying spectra with SNR near 2, using exposure times as low as 0.001 seconds. The methodology involves averaging numerous low-SNR spectra to create augmented databases, then applying unsupervised clustering (autoencoder, PCA, DBSCAN) to filter noise from signals. This approach enables rapid analysis of nano plastics and other low-concentration substances while reducing hardware requirements.

#### D. SERS with Metal-Phenolic Networks and Machine Learning

Surface-enhanced Raman spectroscopy (SERS) combined with metal-phenolic networks (MPNs) and machine learning represents a cutting-edge approach for nano plastic detection. MPNs concentrate nano plastics from low-abundance samples, enabling analysis via full-spectrum SERS. Tailored ML pipelines for outlier filtering, classification, and quantification allow accurate polymer identification and concentration prediction. Detection limits below 0.1 ppm have been achieved, with applications extending to food matrices such as leafy vegetables.




Fig 3: Analysis of SERS

A novel SERS sensor developed through in situ synthesis of gold-core silver-shell nanoparticles on bacterial cellulose (BC@Au@Ag) enables reliable detection of polyethylene and polystyrene MNPs in kale samples at  $4 \text{ mg kg}^{-1}$ , with theoretical detection limits as low as  $1.22 \text{ mg kg}^{-1}$  for PS. The bacterial cellulose provides both a biocompatible scaffold and eco-friendly reducing agent, facilitating green synthesis of nanoparticles with uniform pore structure for reproducible SERS activity.

#### E. Electrochemical and Microfluidic Platforms

Emerging electrochemical sensing strategies, including electrochemical impedance spectroscopy (EIS), photoelectrochemical (PEC), and voltametric sensors, exploit polymer-electrode interfacial interactions for label-free, sensitive, and real-time detection. These platforms offer potential for miniaturization and in-field deployment. Microfluidic-based capture systems and lab-on-chip devices are also advancing toward automated MNP detection in environmental and biological matrices.

#### F. Challenges

Despite these advances, significant challenges remain. Matrix effects from complex environmental and food samples interfere with detection. The lack of certified reference materials for NPs hampers method validation and interlaboratory comparison. Contamination risks during sample preparation can lead to false positives. Distinguishing NPs from natural colloids (e.g., organic matter, clay particles) remains difficult due to similar size and behaviour. Spectral overlap between polymer types and interference from co-contaminants complicate identification. Recent reviews emphasize the need for multimodal workflows and machine learning to improve throughput and specificity.

### V. METHODOLOGY

#### A. Dataset & Model

Combined public MP datasets → 7,200 images (fibre, film, fragment, pellet).  
 YOLOv11m trained at 1280 px, 300 epochs.  
 Post-processing: size calibration (user-provided magnification), statistics export.

#### B. Web Architecture

Frontend: React (drag-and-drop, responsive dashboard).  
 Backend: Fast API + Ultralytics inference.  
 Deployment: Docker on Render/Railway/Hugging Face (free tier viable).

C. Results

Quantitative Performance (500-image test set)

Model	YOLO11m
mAp@50	95.4%
mAp@50:95	71.4%
Precision	94.2%
Recall	92.8%
CPU Inference	51 Ms
GPU Inference	6.4 Ms

Morphology-specific AP@50 >94% across classes. Counting accuracy: 96.8%; size MAE <7% with correct scale.




Fig 4: Workflow

VI. DISCUSSION

A. Strengths of MPWebAI

Zero-cost analysis phase; accessible worldwide.

Sub-10 s response on modest cloud instances.

Immediate quantitative insights for citizen-science and regulatory use and addresses the need for scalable, automated monitoring systems that can operate beyond specialized laboratories.

B. Limitations for Microplastics

Reduced performance on very dense or low-contrast fields. Morphology only (no polymer chemistry)—a limitation shared with many imaging-based approaches.

C. Nano plastics-Specific Limitations & Opportunities

Current optical-microscopy + YOLO pipeline cannot reliably detect particles <200–500 nm. Even high-magnification optical images suffer from diffraction limits. However, emerging technologies offer integration pathways:

**Support for advanced imaging inputs:** The platform could accept images from optical sieves, SEM/TEM, or nano-DIHM systems, with retrained YOLOv11-seg models for NP detection.

**Hybrid pipeline with fluorescence/Raman-enhanced images:** Integration of SERS data or fluorescence staining (e.g., Nile red) could provide chemical specificity.

**Optical sieve module:** Low-cost capture and analysis using test strips readable by standard microscopes could be paired with MPWebAI's quantification engine.

**Multi-modal fusion with spectral data:** Machine learning models that combine imaging features with Raman or FTIR spectra could improve classification accuracy and enable polymer identification.

**AI-enhanced spectroscopy integration:** The platform could incorporate the rapid, low-SNR Raman analysis enabled by Bi-CNNs, extending detection capabilities to nanoparticles with minimal hardware requirements.

These extensions position MPWebAI as a unified platform evolving from MPs to full MNPs monitoring, addressing the analytical challenges of nano plastic traceability, spectral overlap, and matrix interference.

## VII. APPLICATIONS

The proposed web-based intelligent microplastic detection system is applicable across environmental monitoring, water quality analysis, and industrial assessment due to its automated and scalable design.

- **Marine and Freshwater Monitoring:** Enables rapid detection of microplastics in oceans, rivers, and lakes using microscopic images, supporting pollution hotspot identification and long-term ecological studies while reducing manual effort. Recent advances in optical sieve technology demonstrate the feasibility of on-site NP analysis in such environments.
- **Drinking Water Quality Assessment:** Assists laboratories in screening bottled, tap, and groundwater samples to ensure regulatory compliance and evaluate filtration effectiveness. Studies have documented widespread MNP contamination in bottled water and beverages, underscoring the need for accessible monitoring tools.
- **Wastewater Treatment Analysis:** Helps treatment plants measure microplastic removal efficiency across processing stages, supporting process optimization and environmental auditing.
- **Environmental and Academic Research:** Supports researchers and students in conducting microplastic distribution studies, morphological analysis, and standardized dataset generation. Open datasets like MiNa support reproducibility and improve development of generalizable AI models.
- **Industrial Monitoring:** Enables plastic manufacturing and recycling industries to monitor microplastic leakage in effluents and maintain environmental compliance.
- **Food Safety Applications:** With appropriate extensions for SERS or spectroscopic analysis, the platform could support detection of MNPs in agricultural products, addressing concerns about food chain contamination.

## VIII. CONCLUSION

MPWebAI delivers a practical, open-source solution for microplastic monitoring today and provides a clear roadmap for nano plastics integration using 2025–2026 breakthroughs in optical sieves, AI-enhanced spectroscopy, and hybrid detection platforms. By lowering barriers to entry, it accelerates global research, policy, and remediation efforts. As nano plastic pollution enters what researchers describe as a new "nano era" of investigation, platforms that democratize access to advanced detection capabilities will be essential for understanding and mitigating the health and environmental risks of these pervasive contaminants.

## REFERENCES

- [1]. Kumar, A., et al. (2026). Recent progress and technological advancements for detection of micro/nano-plastics in the environment. *Advances in Colloid and Interface Science*, 351, 103817.
- [2]. University of Stuttgart. (2025). New and simple detection method for nano plastics. *EurekAlert! / Nature Photonics*.
- [3]. National Research Council Canada. (2025). In-situ and real-time detection of micro/nanoplastics in water: Combining laboratory experiments and modelling studies. *SciProfiles*.
- [4]. Kousheh, S., Mustapha, A., & Lin, M. (2026). In situ synthesis of gold-core silver-shell nanoparticles on bacterial cellulose for SERS detection of micro- and nanoplastic particles in vegetables. *Analyst*, Advance Article.
- [5]. Lim, J., & Shin, D. (2025). Machine learning-enhanced Raman spectroscopy for fast nanoplastic detection at low SNR. *Sensors and Actuators B: Chemical*, 417, 136892.
- [6]. Das, S., Gebremeskal, Y.H., Maksimova, B.V., Zun, P., & Eremeeva, N.B. (2025). Microplastic contamination and detection in food systems: a review of machine learning, traditional methods, and other relevant factors. *Microchemical Journal*, 218, 115440.



- [7]. Ma, C., Cai, H., Su, L., Chen, Q., & Shi, H. (2026). Advancements in environmental nanoplastics research. *Chinese Science Bulletin*, 71(2), 479-493.
- [8]. University of Melbourne. (2025). Nanoplastics detection chip revolutionizes plastic pollution monitoring. *EurekAlert! / Nature Photonics*.
- [9]. Oliveira, M.J.S., Constantino, C.J.L., & Yang, T. (2026). Nanoplastics as Carriers of Synthetic Dyes and Their Detection Using Metal-Phenolic Networks-enabled SERS and Machine Learning. *FAPESP*.
- [10]. Talanta Open. (2025). Emerging analytical frontiers in microplastic detection: From spectroscopy to smart sensor technologies. *Talanta Open*, 12, 100514.
- [11]. Kumar, A., et al. (2025). Deep-learning enabled rapid and low-cost detection of microplastics in consumer products following on-site extraction and image processing. *RSC Advances*.
- [12]. Vengatesh, T., et al. (2025). A Machine Learning Approach To Microplastic Detection And Quantification In Aquatic Environments. *International Journal of Environmental Sciences*.